


Vials Rack


Acrylic vial racks are resistant to most solvents.
Transparent racks easy for cleaning.
Four sizes of racks suitable for all sample vials.



P/N	Size(mm)	Packing	Description
VR02	18*90*25	1	Vial Rack, 50#, Hole Dia. 12.5mm for 2ml Vial
VR04	230*112*30	1	Vial Rack, 50#, Hole Dia. 16.5mm for 4ml Vials
VR20	290*140*35	1	Vial Rack, 50#, Hole Dia. 22.5mm for Headspace Vials
VR24	350*173*45	1	Vial Rack, 50#, Hole Dia. 28.5mm for 24-400 Vials

Crimper and Decrimper

Crimpers provide a secure vial closure.
Manual decrimpers allow easy remove of aluminum seals without breakage.



P/N	Use	Packing	Description
CR11	Attaches 11mm aluminum crimp seals.	1	11mm Manual Crimper, stainless steel
CR20	Attaches 20mm aluminum crimp seals.	1	20mm Manual Crimper, stainless steel
DCR20	Removes 20mm aluminum crimp seals.	1	20mm Manual Decrimper, stainless steel

HPLC Syringes



P/N	Description
SY01	Disposable HPLC Syringe, 1ml, Slip
SY03	Disposable HPLC Syringe, 3ml, Slip
SY05	Disposable HPLC Syringe, 5ml, Slip
SY10	Disposable HPLC Syringe, 10ml, Slip
SY20	Disposable HPLC Syringe, 20ml, Slip
SY01L	Disposable HPLC Syringe, 1ml, Luer Lock
SY03L	Disposable HPLC Syringe, 3ml, Luer Lock
SY05L	Disposable HPLC Syringe, 5ml, Luer Lock
SY10L	Disposable HPLC Syringe, 10ml, Luer Lock
SY20L	Disposable HPLC Syringe, 20ml, Luer Lock